opti-fix

TAILSTOCKS AND ACCESSORIES	APPLICA	ATION PRODUCT LINES
Part No.: K550879 Tailstock with centre Tip height: 20 mm Adjust. range: 16 mm For receiving of specimens between the tips.	Accordance	opti-fix
Part No.: K550333 Tailstock centre stepped for tailstock (K550879) Tip: Ø 6 mm The step at the tip makes length scanning easier, especially in case of the transmitted light method.	36	opti-fix
Part No.: K550986 Tailstock with spring-loaded centre Tip height: 20 mm Stroke: 8 mm In combination with part No. 550879 for easy and quick receiving of the specimens between the tips.	20 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	opti-fix
Part No.: K550334 Tailstock centre stepped for tailstock with spring-loaded centre (K550986) Tip: Ø 10 mm The step at the tip makes length scanning easier, especially in case of the transmitted light method.	\$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$	opti-fix
Part No.: K550994 Centering unit For centering of parts without centre hole.	2 15	opti-fix